

Anti-LINC00310 antibody (STJ196606)

STJ196606

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Putative Uncharacterized Protein Encoded By Linc00310 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	LINC00310
Uniprot ID	CU082_HUMAN
Immunogen	Synthesized peptide derived from human CU082
Immunogen Region	
Specificity	LINC00310 polyclonal antibody (Putative Uncharacterized Protein Encoded By Linc00310) binds to endogenous Putative Uncharacterized Protein Encoded By Linc00310.
Immunogen Sequence	