

## Anti-mCherry-Tag antibody [6B3] (STJ196637)

STJ196637

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Mouse monoclonal antibody anti-mCherry-Tag is suitable for use in Western Blot research applications. |
| <b>Applications</b>      | WB  |
| <b>Host/Source</b>       | Mouse   |
| <b>Reactivity</b>        | Tag   |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Clonality</b>           | Monoclonal   |
| <b>Clone ID</b>            | 6B3  |
| <b>Concentration</b>       |  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Purification</b>        | The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads. |
| <b>Dilution Range</b>      | WB 1:5000  |
| <b>Formulation</b>         | PBS, pH 7.4, 0.5% BSA, 0.02% Sodium Azide and 50% Glycerol.  |
| <b>Isotype</b>             | IgG1   |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                         |

### TARGET INFORMATION

|                           |   |
|---------------------------|---|
| <b>Gene ID</b>            |   |
| <b>Gene Symbol</b>        |   |
| <b>Uniprot ID</b>         |   |
| <b>Immunogen</b>          | Recombinant Protein of mCherry-Tag                |
| <b>Immunogen Region</b>   |   |
| <b>Specificity</b>        | The antibody detects mCherry tag fusion proteins. |
| <b>Immunogen Sequence</b> |   |