

Anti-Phospho-BAD-S112 antibody (STJ22032)

STJ22032

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Bad-S112 is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	572
Gene Symbol	BAD
Uniprot ID	BAD_HUMAN
Immunogen	A phospho specific peptide corresponding to residues surrounding S112 of human Bad
Immunogen Region	
Specificity	
Immunogen Sequence	