

Anti-CTSG antibody (21-255) (STJ23282)

STJ23282

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CTSG (21-255) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1511
Gene Symbol	CTSG
Uniprot ID	CATG_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-255 of human CTSG (NP_001902.1).
Immunogen Region	21-255
Specificity	
Immunogen Sequence	