

Anti-FUT6 antibody (170-359) (STJ23720)

STJ23720

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FUT6 (170-359) is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:100-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2528
Gene Symbol	FUT6
Uniprot ID	FUT6_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-359 of human FUT6 (NP_000141.1).
Immunogen Region	170-359
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin-embedded human liver using FUT6 antibody (STJ23720) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human stomach using FUT6 antibody (STJ23720) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081