

## Anti-Tri-Methyl-Histone H3-K9 antibody (STJ23986)

STJ23986

### GENERAL INFORMATION


|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Tri-Methyl-Histone H3-K9 is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence. |
| <b>Applications</b>      | WB, IHC, IF  |
| <b>Host/Source</b>       | Rabbit   |
| <b>Reactivity</b>        | Human, Mouse, Rat  |

### PRODUCT PROPERTIES


|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       |   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Purification</b>        | Affinity purification   |
| <b>Dilution Range</b>      | DB 1:500-1:2000<br>WB 1:500-1:2000<br>IHC 1:50-1:200<br>IF 1:50-1:200 |
| <b>Formulation</b>         | PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.               |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store in a freezer at -20°C and avoid freeze-thaw cycles.             |

### TARGET INFORMATION


|                           |   |
|---------------------------|---|
| <b>Gene ID</b>            | 8290  |
| <b>Gene Symbol</b>        | H3-4  |
| <b>Uniprot ID</b>         | H31T_HUMAN  |
| <b>Immunogen</b>          | A synthetic methylated peptide corresponding to residues surrounding K9 of human histone H3 |
| <b>Immunogen Region</b>   |   |
| <b>Specificity</b>        |   |
| <b>Immunogen Sequence</b> |   |




Western blot analysis of extracts of various cell lines, using TriMethyl-Histone H3-K9 antibody (STJ23986) at 1:500 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Dot-blot analysis of all sorts of methylation peptides using TriMethyl-Histone H3-K9 antibody (STJ23986).



Immunohistochemistry of paraffin-embedded rat testis using TriMethyl-Histone H3-K9 antibody (STJ23986) at dilution of 1:200 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry of paraffin-embedded mouse testis using TriMethyl-Histone H3-K9 antibody (STJ23986) at dilution of 1:200 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081