

Anti-MAP3K5 antibody (1173-1374) (STJ24478) STJ24478

GENERAL INFORMATION

Product Type Primary antibodies Short Description Rabbit polyclonal antibody anti-ASK1 (1173-1374) is suitable for use in Immunohistochemistry. Applications IHC Host/Source Rabbit Reactivity Human, Rat

PRODUCT PROPERTIES

Clonality Polyclonal Clone ID Concentration Conjugation Unconjugated Purification Affinity purification Dilution Range IHC 1:50-1:200 Formulation PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3. Isotype IgG Storage Instruction Store in a freezer at-20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID 4217 Gene Symbol MAP3K5 Uniprot ID M3K5_HUMAN (NP_005914.1). Immunogen Region 1173-1374 Specificity Immunogen

Sequence

Immunogen Recombinant fusion protein containing a sequence corresponding to amino acids 1173-1374 of human ASK1

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes. St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081