

Anti-MCAM antibody (447-646) (STJ24519)

STJ24519

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MCAM (447-646) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4162
Gene Symbol	MCAM
Uniprot ID	MUC18_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 447-646 of human CD146/CD146/MCAM (NP_006491.2).
Immunogen Region	447-646
Specificity	
Immunogen Sequence	



Western blot analysis of extracts of HeLa cells, using CD146/CD146/MCAM antibody (STJ24519) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.