

Anti-PRDX4 antibody (38-271) (STJ25106)

STJ25106

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PRDX4 (38-271) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10549
Gene Symbol	PRDX4
Uniprot ID	PRDX4_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 38-271 of human Peroxiredoxin 4 (Peroxiredoxin 4 (PRDX4)) (NP_006397.1).
Immunogen Region	38-271
Specificity	
Immunogen Sequence	



Western blot analysis of extracts of various cell lines, using Peroxiredoxin 4 (Peroxiredoxin 4 (PRDX4)) antibody (STJ25106) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nontat dry milk in TBST.



Western blot analysis of extracts from normal (control) and Peroxiredoxin 4 (Peroxiredoxin 4 (PRDX4)) knockout (KO) 293T cells, using Peroxiredoxin 4 (Peroxiredoxin 4 (PRDX4)) antibody (STJ25106) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nontat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081