

Anti-RAP2B antibody (1-183) (STJ25300)

STJ25300

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RAP2B (1-183) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:100-1:500
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5912
Gene Symbol	RAP2B
Uniprot ID	RAP2B_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-183 of human RAP2B (NP_002877.2).
Immunogen Region	1-183
Specificity	
Immunogen Sequence	



Western blot analysis of extracts of various cell lines, using RAP2B antibody (STJ25300) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.