

Anti-RETN antibody (1-108) (STJ25340)

STJ25340

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RETN (1-108) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56729
Gene Symbol	RETN
Uniprot ID	RETN_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-108 of human RETN (NP_001180303.1).
Immunogen Region	1-108
Specificity	
Immunogen Sequence	