

Anti-TCTN3 antibody (200-475) (STJ25800)

STJ25800

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TCTN3 (200-475) is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26123
Gene Symbol	TCTN3
Uniprot ID	TECT3_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-475 of human TCTN3 (NP_056446.4).
Immunogen Region	200-475
Specificity	
Immunogen Sequence	



Western blot analysis of extracts of various cell lines, using TCTN3 antibody (STJ25800) at 1:500 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST.