

Anti-TFF3 antibody (22-94) (STJ26190)

STJ26190

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TFF3 (22-94) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7033
Gene Symbol	TFF3
Uniprot ID	TFF3_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-94 of human TFF3 (NP_003217.3).
Immunogen Region	22-94
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of human colon carcinoma using TFF3 rabbit polyclonal antibody (STJ26190) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081