

Anti-AZGP1 antibody (21-298) (STJ27318)

STJ27318

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-AZGP1 (21-298) is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	563
Gene Symbol	AZGP1
Uniprot ID	ZA2G_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-298 of human Zinc-Alpha 2-glycoprotein (ZAG/Zinc-Alpha 2-glycoprotein (ZAG/AZGP1)) (NP_001176.1).
Immunogen Region	21-298
Specificity	
Immunogen Sequence	



Western blot analysis of extracts of human serum, using Zinc-Alpha 2-glycoprotein (ZAG/Zinc-Alpha 2-glycoprotein (ZAG/AZGP1)) antibody (STJ27318) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081