

Anti-Phospho-BTK-Y223 antibody (STJ27514)

STJ27514

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Btk-Y223 is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:20-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	695
Gene Symbol	BTK
Uniprot ID	BTK_HUMAN
Immunogen	A synthetic phosphorylated peptide around Y223 of human Btk (NP_000052.1).
Immunogen Region	
Specificity	
Immunogen Sequence	