

Anti-Phospho-ESPL1-S1126 antibody (STJ27522)

STJ27522

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-ESPL1-S1126 is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:20-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9700
Gene Symbol	ESPL1
Uniprot ID	ESPL1_HUMAN
Immunogen	A synthetic phosphorylated peptide around S1126 of human ESPL1 (NP_036423.4).
Immunogen Region	
Specificity	
Immunogen Sequence	