

Anti-RABL2A antibody (STJ27545)

STJ27545

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Guinea Pig polyclonal antibody anti-RABL2A is suitable for use in Western Blot.
Applications	WB
Host/Source	Guinea pig
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	68708
Gene Symbol	Rabl2
Uniprot ID	RBL2A_MOUSE
Immunogen	Recombinant Protein of mouse RABL2A
Immunogen Region	
Specificity	
Immunogen Sequence	