

Anti-TNFRSF13B antibody (184-293) (STJ27625)

STJ27625

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TNFRSF13B (184-293) is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23495
Gene Symbol	TNFRSF13B
Uniprot ID	TR13B_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 184-293 of human TNFRSF13B (NP_036584.1).
Immunogen Region	184-293
Specificity	
Immunogen Sequence	