

Anti-MRPL28 antibody (1-256) (STJ27702)

STJ27702

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MRPL28 (1-256) is suitable for use in Western Blot, Immunohistochemistry, Immunofluorescence and Immunoprecipitation.
Applications	WB, IHC, IF, IP
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:200 IF 1:50-1:100 IP 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage	Store in a freezer at -20°C and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	10573
Gene Symbol	MRPL28
Uniprot ID	RM28_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-256 of human MRPL28 (NP_006419.2).
Immunogen Region	1-256
Specificity	
Immunogen Sequence	




Western blot analysis of extracts of various cell lines, using MRPL28 antibody (STJ27702) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit.



Immunohistochemistry of paraffin-embedded rat ovary using MRPL28 antibody (STJ27702) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human lung cancer using MRPL28 antibody (STJ27702) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human gastric cancer using MRPL28 antibody (STJ27702) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081