

Anti-PRPF4 antibody (1-260) (STJ27846)

STJ27846

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PRPF4 (1-260) is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION


Gene ID	9128
Gene Symbol	PRPF4
Uniprot ID	PRP4_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human PRPF4 (NP_004688.2).
Immunogen Region	1-260
Specificity	
Immunogen Sequence	




Western blot analysis of extracts of various cell lines, using PRPF4 antibody (STJ27846) at 1:1000 dilution. Secondary antibody: HRP-Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.



Immunohistochemistry of paraffin-embedded rat lung using PRPF4 rabbit polyclonal antibody (STJ27846) at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded rat spleen using PRPF4 rabbit polyclonal antibody (STJ27846) at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded human brain using PRPF4 rabbit polyclonal antibody (STJ27846) at dilution of 1:200 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081