

Anti-TFPI2 antibody (23-235) (STJ27970)

STJ27970

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TFPI2 (23-235) is suitable for use in Western Blot and Immunohistochemistry.
Applications	WB, IHC
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7980
Gene Symbol	TFPI2
Uniprot ID	TFPI2_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-235 of human TFPI2 (NP_006519.1).
Immunogen Region	23-235
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin-embedded human liver cancer using TFPI2 antibody (STJ27970) at dilution of 1:200 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081