

Anti-BAG5 antibody (198-447) (STJ28358)

GENERAL INFORMATION

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-BAG5 (198-447) is suitable for use in Western Blot and Immunohistochemistry.

Applications WB, IHC Host/Source Rabbit

Reactivity Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality Polyclonal

Clone ID

Concentration

Conjugation Unconjugated Purification Affinity purification Dilution Range WB 1:500-1:2000 IHC 1:50-1:200

Formulation PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.

Isotype IgG

Storage Instruction Store in a freezer at-20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID 9529 Gene Symbol BAG5

Uniprot ID BAG5_HUMAN

Immunogen Recombinant fusion protein containing a sequence corresponding to amino acids 198-447 of human BAG5 (NP_004864.1).

Immunogen Region 198-447 Specificity Immunogen Sequence





