

Anti-XCL1 antibody (22-114) (STJ28491)

STJ28491

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-XCL1 (22-114) is suitable for use in Immunofluorescence. |
| Applications | IF |
| Host/Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | IF 1:50-1:100 |
| Formulation | PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3. |
| Isotype | IgG |
| Storage Instruction | Store in a freezer at -20°C and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 6375 |
| Gene Symbol | XCL1 |
| Uniprot ID | XCL1_HUMAN |
| Immunogen | Recombinant fusion protein containing a sequence corresponding to amino acids 22-114 of human XCL1 (NP_002986.1). |
| Immunogen Region | 22-114 |
| Specificity | |
| Immunogen Sequence | |