

Anti-C1D antibody (1-141) (STJ28531)

STJ28531

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-C1D (1-141) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10438
Gene Symbol	C1D
Uniprot ID	C1D_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-141 of human C1D (NP_001177194.1).
Immunogen Region	1-141
Specificity	
Immunogen Sequence	