

Anti-TFPT antibody (1-253) (STJ28544)

STJ28544

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TFPT (1-253) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29844
Gene Symbol	TFPT
Uniprot ID	TFPT_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-253 of human TFPT (NP_037474.1).
Immunogen Region	1-253
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of U2OS cells using TFPT antibody (STJ28544).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081