

Anti-HPS1 antibody (1-324) (STJ28704)

STJ28704

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-HPS1 (1-324) is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence.
Applications	WB, IHC, IF
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:100 IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3257
Gene Symbol	HPS1
Uniprot ID	HPS1_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-324 of human HPS1 (NP_872577.1).
Immunogen Region	1-324
Specificity	
Immunogen Sequence	