

Anti-RFX5 antibody (367-616) (STJ29048)

STJ29048

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RFX5 (367-616) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5993
Gene Symbol	RFX5
Uniprot ID	RFX5_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 367-616 of human RFX5 (NP_000440.1).
Immunogen Region	367-616
Specificity	
Immunogen Sequence	