

Anti-CAMK1 antibody (271-370) (STJ29098)

STJ29098

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CAMK1 (271-370) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8536
Gene Symbol	CAMK1
Uniprot ID	KCC1A_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 271-370 of human CaMKI (NP_003647.1).
Immunogen Region	271-370
Specificity	
Immunogen Sequence	



Immunofluorescence analysis of U2OS cells using CaMK1 antibody (STJ29098). Blue: DAPI for nuclear staining.