

Anti-NAT9 antibody (1-207) (STJ29284)

STJ29284

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NAT9 (1-207) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26151
Gene Symbol	NAT9
Uniprot ID	NAT9_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-207 of human NAT9 (NP_056469.2).
Immunogen Region	1-207
Specificity	
Immunogen Sequence	