

Anti-ZDHHC16 antibody (207-296) (STJ29306)

STJ29306

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ZDHHC16 (207-296) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84287
Gene Symbol	ZDHHC16
Uniprot ID	ZDH16_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 207-296 of human ZDHHC16 (NP_932162.1).
Immunogen Region	207-296
Specificity	
Immunogen Sequence	