

Anti-Phospho-Histone H3-S32 antibody (STJ29330)

STJ29330

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Histone H3.3-S32 is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3020/3021
Gene Symbol	H3-3A.H3-3B
Uniprot ID	H33_HUMAN
Immunogen	A phospho specific peptide corresponding to residues surrounding S32 of human Histone H3.3
Immunogen Region	
Specificity	
Immunogen Sequence	