

Anti-Phospho-KSR1-S392 antibody (STJ29331)

STJ29331

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-KSR1-S392 is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8844
Gene Symbol	KSR1
Uniprot ID	KSR1_HUMAN
Immunogen	A phospho specific peptide corresponding to residues surrounding S392 of human KSR1
Immunogen Region	
Specificity	
Immunogen Sequence	