

Anti-DPF1 antibody (295-414) (STJ29488)

STJ29488

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-DPF1 (295-414) is suitable for use in Western Blot and Immunofluorescence. |
| Applications | WB, IF |
| Host/Source | Rabbit |
| Reactivity | Human, Mouse, Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | WB 1:500-1:2000 IF 1:50-1:100 |
| Formulation | PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3. |
| Isotype | IgG |
| Storage Instruction | Store in a freezer at -20°C and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 8193 |
| Gene Symbol | DPF1 |
| Uniprot ID | DPF1_HUMAN |
| Immunogen | Recombinant fusion protein containing a sequence corresponding to amino acids 295-414 of human DPF1 (NP_001128627.1). |
| Immunogen Region | 295-414 |
| Specificity | |
| Immunogen Sequence | |