

Anti-NAA40 antibody (1-237) (STJ29522)

STJ29522

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NAA40 (1-237) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79829
Gene Symbol	NAA40
Uniprot ID	NAA40_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-237 of human NAA40 (NP_079047.2).
Immunogen Region	1-237
Specificity	
Immunogen Sequence	