

Anti-PRKCD antibody (477-676) (STJ29804)

STJ29804

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PRKCD (477-676) is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence.
Applications	WB, IHC, IF
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IHC 1:50-1:100 IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5580
Gene Symbol	PRKCD
Uniprot ID	KPCD_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 477-676 of human PKC delta (NP_997704.1).
Immunogen Region	477-676
Specificity	
Immunogen Sequence	




Western blot analysis of extracts of various cell lines, using PKC Delta antibody (STJ29804) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1min.



Immunohistochemistry of paraffin-embedded rat lung using PKC Delta antibody (STJ29804) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human tonsil using PKC Delta antibody (STJ29804) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse testis using PKC Delta antibody (STJ29804) at dilution of 1:100 (40x lens).