

Anti-AFP antibody [5108] (STJ400006)

STJ400006

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-AFP is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5108
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Purity accurate to greater than or equal to 95%
Dilution Range	
Formulation	0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Isotype	IgG1
Storage Instruction	24 months from manufacturing at 2°C. Temperature analysis:-70 °C, 21 days N/D-20 °C, 21 days N/D +4 °C, 21 days N/D +35 °C, 7 days N/D +35 °C, 21 days N/D +45 °C, 3 days N/D +45 °C, 7 days N/D. Shelf life stability testing at 2-8 °C in the product bu

TARGET INFORMATION

Gene ID	174
Gene Symbol	AFP
Uniprot ID	FETA_HUMAN
Immunogen	human Alpha1-fetoprotein
Immunogen Region	
Immunogen Specificity	
Immunogen Sequence	