

Anti-Alpha Subunit antibody [6601] (STJ400011)

STJ400011

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Alpha Subunit is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	6601
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Purity accurate to greater than or equal to 95%
Dilution Range	
Formulation	37 mM citrate, 125 mM phosphate, pH 6.0, 0.095 % NaN ₃ as a preservative
Isotype	IgG2a
Storage	36 months from manufacturing at 2°C. Temperature analysis:-70 °C, 21 days N/D-20 °C, 21 days OK +4 °C, 21 days OK +30 °C, 21 days OK +35 °C, 21 days OK +45 °C, 7 days OK. Shelf life stability testing at 2-8 °C in the product buffer.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	the alpha subunit of human chorionic gonadotropin (hCG) , luteinizing hormone (LH) , follicle-stimulating hormone (FSH) , and thyroid-stimulating hormone (TSH)
Immunogen Region	
Specificity	
Immunogen Sequence	