

Anti-CA125 antibody [4601] (STJ400019)

STJ400019

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CA125 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4601
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Purity accurate to greater than or equal to 95%
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage	Unspecified, storage at 2°C. Temperature analysis:-70 °C, 21 days OK-20 °C, 21 days OK +4 °C, 21 days OK +35 °C, 7 days OK +35 °C, 21 days OK +45 °C, 3 days OK +45 °C, 7 days OK. Shelf life stability testing at 2–8 °C in the product buffer.
Instruction	

TARGET INFORMATION

Gene ID	94025
Gene Symbol	MUC16
Uniprot ID	MUC16_HUMAN
Immunogen	human cancer antigen 125 (CA125)
Immunogen Region	
Specificity	
Immunogen Sequence	