

Anti-AVP antibody [4806] (STJ400044)

STJ400044

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-AVP is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4806
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Purity accurate to greater than or equal to 95%
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Unspecified, storage at 2°C. Temperature analysis: -70 °C, 21 days OK -20 °C, 21 days OK +4 °C, 21 days OK +35 °C, 7 days OK +35 °C, 21 days OK +45 °C, 3 days OK +45 °C, 7 days OK. Shelf life stability testing at 2-8 °C in the product buffer.

TARGET INFORMATION

Gene ID	551
Gene Symbol	AVP
Uniprot ID	NEU2_HUMAN
Immunogen	human copeptin
Immunogen Region	
Specificity	
Immunogen Sequence	