

Anti-TNNI3 antibody [9705] (STJ400053)

STJ400053

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TNNI3 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	9705
Concentration	1.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Purity accurate to greater than or equal to 95%
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Isotype	IgG1
Storage	36 months from manufacturing at 2°C. Temperature analysis:-70 °C, 21 days N/D-20 °C, 21 days OK +4 °C, 21 days OK +30 °C, 21 days OK +35 °C, 21 days OK +45 °C, 7 days OK. Shelf life stability testing at 2-8 °C in the product buffer.
Instruction	

TARGET INFORMATION

Gene ID	7137
Gene Symbol	TNNI3
Uniprot ID	TNNI3_HUMAN
Immunogen	human cardiac troponin I (cTnI)
Immunogen Region	
Immunogen Specificity	
Immunogen Sequence	