

## Anti-GH1 antibody [5801] (STJ400071)

STJ400071

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-GH1 is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	5801
<b>Concentration</b>	1.0 mg/ml (+/-10 %)
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Purity accurate to greater than or equal to 95%
<b>Dilution Range</b>	
<b>Formulation</b>	0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	36 months from manufacturing at 2°C. Temperature analysis:-70 °C, 21 days N/D-20 °C, 21 days N/D +4 °C, 21 days N/D +35 °C, 7 days N/D +35 °C, 21 days N/D +45 °C, 3 days N/D +45 °C, 7 days N/D. Shelf life stability testing at 2-8 °C in the product bu

### TARGET INFORMATION

<b>Gene ID</b>	2688
<b>Gene Symbol</b>	GH1
<b>Uniprot ID</b>	SOMA_HUMAN
<b>Immunogen</b>	human growth hormone (GH)
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	