

Anti-CG beta antibody [5008] (STJ400083)

STJ400083

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CG beta is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5008
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Purity accurate to greater than or equal to 95%
Dilution Range	
Formulation	0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Isotype	IgG1
Storage Instruction	36 months from manufacturing at 2°C

TARGET INFORMATION

Gene ID	1082/93659/94115
Gene Symbol	CGB3.CGB5.CGB8
Uniprot ID	CGB3_HUMAN
Immunogen	human chorionic gonadotropin and its free beta subunit
Immunogen Region	
Specificity	
Immunogen Sequence	