

## Anti-CG antibody [5016] (STJ400087)

STJ400087

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CG is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	5016
<b>Concentration</b>	5.0 mg/ml (+/-10 %)
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Purity accurate to greater than or equal to 95%
<b>Dilution Range</b>	
<b>Formulation</b>	37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	24 months from manufacturing at 2°C

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">1511</a>
<b>Gene Symbol</b>	<a href="#">CTSG</a>
<b>Uniprot ID</b>	<a href="#">CATG_HUMAN</a>
<b>Immunogen</b>	human chorionic gonadotropin and its free beta subunit
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	