

Anti-LCN2 antibody [4204] (STJ400140)

STJ400140

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-LCN2 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4204
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Purity accurate to greater than or equal to 95%
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage Instruction	Unspecified, storage at 2°C. Temperature analysis: -70 °C, 21 days N/D-20 °C, 21 days OK +4 °C, 21 days OK +30 °C, 21 days OK +35 °C, 7 days OK +35 °C, 21 days Failed due to reduced antigen binding +45 °C, 3 days OK +45 °C, 7 days OK, some charge alte

TARGET INFORMATION

Gene ID	3934
Gene Symbol	LCN2
Uniprot ID	NGAL_HUMAN
Immunogen	human neutrophil gelatinase-associated lipocalin (NGAL)
Immunogen Region	
Immunogen Specificity	
Immunogen Sequence	