

Anti-TSHB antibody [5408] (STJ400196)

STJ400196

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TSHB is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5408
Concentration	1.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Purity accurate to greater than or equal to 95%
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage	36 months from manufacturing at 2°C. Temperature analysis:-70 °C, 21 days N/D-20 °C, 21 days N/D +4 °C, 21 days N/D +35 °C, 7
Instruction	days N/D +35 °C, 21 days N/D +45 °C, 3 days N/D +45 °C, 7 days N/D. Shelf life stability testing at 2-8 °C in the product bu

TARGET INFORMATION

Gene ID	7252
Gene Symbol	TSHB
Uniprot ID	TSHB_HUMAN
Immunogen	human thyrotropin and its free beta subunit
Immunogen Region	
Specificity	
Immunogen Sequence	