

## Anti-NT-proBNP antibody (22-34) [1308] (STJ400230)

STJ400230

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-NT-proBNP (22-34) is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	1308
<b>Concentration</b>	5.0 mg/ml (+/-10 %)
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Greater than 95% purity
<b>Dilution Range</b>	
<b>Formulation</b>	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
<b>Isotype</b>	IgG2b
<b>Storage Instruction</b>	Unspecified, storage at 2–8°C. Temperature analysis: -70 °C, 21 days Not Determined (N/D)-20 °C, 14 days OK +4 °C, 14 days OK +35 °C, 14 days OK +45 °C, 7 days OK. Shelf life stability testing at 2–8 °C in the product buffer.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Human N-terminal pro-B-type natriuretic peptide
<b>Immunogen Region</b>	22-34
<b>Specificity</b>	
<b>Immunogen Sequence</b>	