

Anti-Galectin-3 antibody (A66-T104) [10303] (STJ400240)

STJ400240

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Galectin-3 (A66-T104) is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	10303
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Greater than 95% purity
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage	Unspecified, storage at 2-8°C. Temperature analysis: -70 °C, 21 days Not Determined (N/D)-20 °C, 14 days OK +4 °C, 14 days OK +35 °C, 14 days OK +45 °C, 7 days OK. Shelf life stability testing at 2-8 °C in the product buffer.
Instruction	

TARGET INFORMATION

Gene ID	3958
Gene Symbol	LGALS3
Uniprot ID	LEG3_HUMAN
Immunogen	Human galectin-3
Immunogen Region	A66-T104
Specificity	
Immunogen Sequence	