

Anti-Lp-PLA2 antibody [11207] (STJ400248)

STJ400248

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Lp-PLA2 is suitable for use in ELISA research applications. |
| Applications | ELISA |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | 11207 |
| Concentration | 5.0 mg/ml (+/-10 %) |
| Conjugation | Unconjugated |
| Purification | Greater than 95% purity |
| Dilution Range | |
| Formulation | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative |
| Isotype | IgG1 |
| Storage Instruction | Unspecified, storage at 2-8°C. Temperature analysis: -70 °C, 21 days Not Determined (N/D)-20 °C, 14 days OK +4 °C, 14 days OK +35 °C, 14 days OK +45 °C, 7 days OK. Shelf life stability testing at 2-8 °C in the product buffer. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | Human lipoprotein-associated phospholipase A2 |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |