

## Anti-AMH antibody (R26-R451) [11302] (STJ400251)

STJ400251

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-AMH (R26-R451) is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	11302
<b>Concentration</b>	5.0 mg/ml (+/-10 %)
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Greater than 95% purity
<b>Dilution Range</b>	
<b>Formulation</b>	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
<b>Isotype</b>	IgG1
<b>Storage</b>	Unspecified, storage at 2-8°C. Temperature analysis: -70 °C, 21 days Not Determined (N/D)-20 °C, 14 days OK +4 °C, 14 days OK +35 °C, 14 days OK +45 °C, 7 days OK. Shelf life stability testing at 2-8 °C in the product buffer.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">268</a>
<b>Gene Symbol</b>	<a href="#">AMH</a>
<b>Uniprot ID</b>	<a href="#">MIS_HUMAN</a>
<b>Immunogen</b>	Human anti-Müllerian hormone (AMH)
<b>Immunogen Region</b>	R26-R451
<b>Specificity</b>	
<b>Immunogen Sequence</b>	