

Anti-AMH antibody (R26-R451) [11303] (STJ400252)

STJ400252

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-AMH (R26-R451) is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	11303
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Greater than 95% purity
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Isotype	IgG1
Storage	12 months from manufacturing at 2–8°C. Temperature analysis:-70 °C, 21 days Not Determined (N/D)-20 °C, 14 days OK +4 °C, 14 days OK +35 °C, 14 days OK +45 °C, 7 days OK. Shelf life stability testing at 2–8 °C in the product buffer.
Instruction	

TARGET INFORMATION

Gene ID	268
Gene Symbol	AMH
Uniprot ID	MIS_HUMAN
Immunogen	Human anti-Müllerian hormone (AMH)
Immunogen Region	R26-R451
Specificity	
Immunogen Sequence	